



# Disposable Medical Mask



*3 Layers Filter Dust* , Pollen, Virus and Bacteria

Middle layer-High Density Nonwoven


**BFE UP 99%**

**50** PCS





# Disposable Medical Mask



*3 Layers Filter Dust* , Pollen, Virus and Bacteria

Middle layer-High Density Nonwoven

**BFE UP 99%**

**50**<sub>pcs</sub>





*Filter Dust, Pollen, Virus and Bacteria*

# Disposable Medical Mask

## 50 pcs

LOT: WK20200420

MFR DATE: 20/04/2020

EXP DATE: 20/04/2022

QUANTITY: 50pcs/box

VENDOR NAME: HuBei Weikang Protective Products Co.Ltd

ADD: GongQing ChuangYe Industrial Zone, YeWang Rd.  
Xiantao 433000, Hubei, China

F&D Administration Registration Number 3005090771





Well Klean

#### Point

- Nonwoven antibiosis material
- Filter by 3 layers
- BFE up 99%

#### Purpose

- Filter harmful bacteria which is bad for humans.

#### Use Method

1



1, The colored is outside, put nose metal strip up.

2



2, Press metal strip on the nose, let mask close to face.

3



3, Wear loops on the ears.

#### Notice for facemask wearing

- Facemask only reduce the risk for virus infections, not guarantee 100% safety.
- Don't adjust the facemask after you finish wearing, and don't wear again when take off the mask. Avoid virus infection from your hands. After take off the mask, should wash your hands.
- When you take off the mask, please don't touch the outside layer of mask, maybe it has virus.
- After take off the mask, put the facemask in the polybag or paperbag, and then put them in the trash can.
- Change the facemask even one day one piece, if facemask broken or dirty, please change it soon.



ISO 13485:2016





*Filter Dust, Pollen, Virus and Bacteria*


# Disposable Medical Mask

## 50 pcs

*3 Layers*

*Earloop*

Middle layer-High  
Density Nonwoven  
BFE UP 99%



Nonwoven